

Ebola virus EBOV (Subtype Sudan, strain Gulu) Glycoprotein / GP1 (mucin domain deleted) Antibody, Rabbit MAb

Catalog Number: 40094-R301



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Sudan ebolavirus EBOV-G / Glycoprotein Protein (Catalog#40094-V08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Sudan ebolavirus EBOV-G / Glycoprotein (rSudan ebolavirus EBOV-G / Glycoprotein; Catalog#40094-V08H; Q7T9D9; Met1-Asp320).
Ig Type:	Rabbit IgG
Clone ID:	301
Specificity:	EBOV-G / Glycoprotein
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	Glycoprotein, GP
APPLICATIONS	
Applications:	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect EBOV-G / Glycoprotein.

Please Note: Optimal concentrations/dilutions should be determined by the end user.