

SARS-COV Spike Antibody, Rabbit MAb

Catalog Number: 40150-R327




Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen	Recombinant SARS-CoV Spike Protein (Catalog#40150-V08B2)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant SARS-CoV Spike (Catalog#40150-V08B2; AAX16192.1; Arg306-Phe527).
Ig Type:	Rabbit IgG
Clone ID:	327
Specificity:	SARS-CoV Spike
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB, ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000

Please Note: Optimal concentrations/dilutions should be determined by the end user.



Anti-SARS-CoV Spike rabbit monoclonal antibody at 1:1000 dilution.

Sample: Recombinant Protein 10 ng

Lane A: SARS-CoV Spike/RBD Protein (RBD, His Tag) (Catalog# 40150-V08B2)

Lane B: SARS-CoV-2 (2019-nCoV) Spike RBD-mFc Recombinant Protein (Catalog# 40592-V05H)

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Manufactured By Sino Biological Inc.

Tel: 215-583-7898 (US); +86-400-890-9989(Global) ●

<http://www.sinobiological.com>