

Influenza A H5N1 (A/chicken/VietNam/NCVD-016/2008) Hemagglutinin (HA1 Subunit) Antibody, Rabbit PAb

Catalog Number: 40158-RP01



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION	
Immunogen:	Recombinant Influenza A H5N1 (A/chicken/VietNam/NCVD-016/2008) Hemagglutinin (HA1 Subunit) Protein (Catalog#40158-V08H1)
Preparation	Produced in rabbits immunized with purified, recombinant Influenza A H5N1 (A/chicken/VietNam/NCVD-016/2008) Hemagglutinin (HA1 Subunit) (Catalog#40158-V08H1; ACO07033.1; Met1-Arg343). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Influenza A H5N1 (A/chicken/VietNam/NCVD-016/2008) Hemagglutinin (HA1 Subunit)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	Hemagglutinin,HA1
APPLICATIONS	
Applications:	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Influenza A H5N1 (A/chicken/VietNam/NCVD-016/2008) Hemagglutinin (HA1 Subunit).

Please Note: Optimal concentrations/dilutions should be determined by the end user.