

HIV-1 p24 Protein (group O, strain BCF06) Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 40256-RP01

GENERAL INFORMATION

Immunogen:	Recombinant HIV-1 p24 Protein (group O, strain BCF06) protein (Catalog#40256-V08E)
Preparation	Produced in rabbits immunized with purified, recombinant HIV-1 p24 Protein (group O, strain BCF06) (Catalog#40256-V08E; BAH97658.1; Pro131-Leu362). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	HIV-1 p24 Protein (group O, strain BCF06)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	Gag-p24

APPLICATIONS

Applications:	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect HIV-1 p24 Protein (group O, strain BCF06).
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.