

Dengue virus DENV-2 (strain New Guinea C) NS5A / Nonstructural protein 5A Antibody, Rabbit PAb

Catalog Number: 40267-RP01



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant DENV Dengue virus DENV-2 (strain New Guinea C) NS5A / Nonstructural protein 5A Protein (Catalog#40267-V07E1)
Preparation	Produced in rabbits immunized with purified, recombinant DENV Dengue virus DENV-2 (strain New Guinea C) NS5A / Nonstructural protein 5A (Catalog#40267-V07E1; AAC59275.1; Gly2492-Gln2787). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	DENV Dengue virus DENV-2 (strain New Guinea C) NS5A / Nonstructural protein 5A
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	NS5A, Nonstructural protein 5A

APPLICATIONS

Applications:	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect DENV Dengue virus DENV-2 (strain New Guinea C) NS5A / Nonstructural protein 5A.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.