

**Ebola virus EBOV (subtype Zaire, strain
H.sapiens-wt/GIN/2014/Kissidougou-C15)
VP24 Antibody, Rabbit PAb, Antigen Affinity
Purified**
Catalog Number: 40454-T46



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) VP24 Protein (Catalog#40454-V07E)
Preparation	Produced in rabbits immunized with purified, recombinant Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) VP24 (Catalog#40454-V07E; AHX24653.1; Met1-Gly233). Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) VP24 specific IgG was purified by Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) VP24 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) VP24
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	VP24

APPLICATIONS

Applications:	WB
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION


Western Blot	WB: 1:1000-1:5000
---------------------	-------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.

**Ebola virus EBOV (subtype Zaire, strain
H.sapiens-wt/GIN/2014/Kissidougou-C15)
VP24 Antibody, Rabbit PAb, Antigen Affinity
Purified**
Catalog Number: 40454-T46



Sino Biological
Biological Solution Specialist



Anti-Ebola virus EBOV (subtype Zaire, strain
H.sapiens-wt/GIN/2014/Kissidougou-C15)
VP24 rabbit polyclonal antibody at 1:2000
dilution.

Sample: Ebola virus EBOV (subtype Zaire,
strain H.sapiens-wt/GIN/2014/Kissidougou-
C15) VP24 Protein

Lane A: 20ng

Lane B: 10ng

Secondary

Goat Anti-Rabbit IgG H&L (Dylight800) at
1/10000 dilution

Performed under reducing conditions.