

AcmNPV Envelope glycoprotein gp64 / AcmNPV-gp64 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 40496-MM02

GENERAL INFORMATION	
Immunogen:	Recombinant AcmNPV Envelope glycoprotein gp64 / AcmNPV-gp64 Protein (Catalog#40496-V08B)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant AcmNPV AcmNPV GP64 (Catalog#40496-V08B; NP_054158.1; Met1-Thr481). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Clone ID:	02
Specificity:	AcmNPV Envelope glycoprotein gp64 / AcmNPV-gp64
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA,ELISA(Cap)
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000
ELISA(Cap)	ELISA(Cap): 1:250-1:2000 This antibody will detect AcmNPV Envelope glycoprotein gp64 / AcmNPV-gp64 in ELISA pair set (Catalog: # SEKA40496). In a sandwich ELISA, it can be used as capture antibody.

Please Note: Optimal concentrations/dilutions should be determined by the end user.