

SARS-CoV-2 Spike Antibody (Biotin), Rabbit Mab



Sino Biological
Biological Solution Specialist

Catalog Number: 40592-R118-B

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	40592-V05H
Preparation	This product is a recombinant monoclonal antibody expressed from HEK293 cells and then biotinylated.
Ig Type:	Rabbit IgG
Clone ID:	118
Specificity:	SARS-CoV-2 Spike RBD Protein (Can block the binding between RBD and ACE2)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	coronavirus s1,coronavirus s2,coronavirus spike,cov spike,ncov RBD,ncov s1,ncov s2,ncov spike,RBD,S,s1,Spike RBD
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 0.1-1 µg/ml

Please Note: Optimal concentrations/dilutions should be determined by the end user.