

Human metapneumovirus (strain HMPV/Homo sapiens/PER/CFI0466/2010/B) glycoprotein G Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 40798-T62

GENERAL INFORMATION

Immunogen:	Recombinant Human metapneumovirus (strain HMPV/Homo sapiens/PER/CFI0466/2010/B) glycoprotein G Protein (Catalog#40798-V08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human metapneumovirus (strain HMPV/Homo sapiens/PER/CFI0466/2010/B) glycoprotein G (Catalog#40798-V08H; AHV79646.1; Asp52-Ser238). Glycoprotein G specific IgG was purified by Human metapneumovirus (strain HMPV/Homo sapiens/PER/CFI0466/2010/B) glycoprotein G affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human metapneumovirus (strain HMPV/Homo sapiens/PER/CFI0466/2010/B) glycoprotein G
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Human metapneumovirus (strain HMPV/Homo sapiens/PER/CFI0466/2010/B) glycoprotein G Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 40798-T62



Anti-Human metapneumovirus glycoprotein G rabbit polyclonal antibody in 2 µg/mL.

Lane A: Human metapneumovirus (strain HMPV/Homo sapiens/PER/CFI0466/2010/B) glycoprotein G Protein (Cat#40798-V08H)(10ng)

Secondary
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution

Developed using the ECL technique.
Performed under reducing conditions.