

Monkeypox Virus (MPXV) Envelope protein H3 Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 40893-T62



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant Monkeypox Virus (MPXV) Envelope protein H3 (Catalog#40893-V08H1)
Preparation	Produced in rabbits immunized with purified, recombinant Monkeypox Virus (MPXV) Envelope protein H3 (Catalog#40893-V08H1; QNP13688.1; Met1-Pro278). H3L Protein specific IgG was purified by Monkeypox Virus (MPXV) Envelope protein H3 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Monkeypox Virus (MPXV) Envelope protein H3
	Has cross-reactivity in WB and ELISA with Monkeypox Virus (MPXV) (Clade Ib) Envelope protein H3L (Cat# 41040-V08H)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:5000-1:50000
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Monkeypox Virus (MPXV) Envelope protein H3 Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 40893-T62



Sino Biological
Biological Solution Specialist




Anti-MPXV-H3 rabbit polyclonal antibody at 1:5000 dilution.

Lane A: Monkeypox Virus (MPXV) Envelope protein H3 (Cat#40893-V08H1) 5ng

Secondary
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution

Developed using the ECL technique.
Performed under reducing conditions.



Anti-MPXV-H3 rabbit polyclonal antibody at 1:5000 dilution.

Lane A: Monkeypox Virus (MPXV) (Clade Ib) Envelope protein H3L (Cat#41040-V08H) 10ng

Secondary
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution

Developed using the ECL technique.
Performed under reducing conditions.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>