

BACE1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50002-T56

GENERAL INFORMATION	
Immunogen:	50002-M08H
Preparation	Produced in rabbits immunized with purified, recombinant Mouse BACE1 (rM BACE1; Catalog#50002-M08H; NP_035922.4; Met1-Thr457). BACE1 specific IgG was purified by Mouse BACE1 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse BACE1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	ASP2
APPLICATIONS	
Applications:	WB,ELISA,IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:500-1:2000
Western Blot	WB: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

BACE1 Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 50002-T56



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of mouse BACE1 in mouse brain with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).



Anti-BACE1 rabbit polyclonal antibody at 1:500 dilution.

Lane A: Mouse brain tissue lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti- Rabbit IgG H&L (Dylight 800) at 1/10000 dilution.

Developed using the Odyssey technique.

Performed under reducing conditions.

Predicted band size: 56 kDa

Observed band size: 70 kDa

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>