

ASAH2 Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50006-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse ASAH2 protein (Catalog#50006-M07H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse ASAH2 Luminal domain (rM ASAH2; Catalog#50006-M07H; Thr 34-Thr 756; NP_061300.1). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse ASAH2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Mouse ASAH2.

Please Note: Optimal concentrations/dilutions should be determined by the end user.