

Carbonic Anhydrase XII / Car12 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50014-R012

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse Carbonic Anhydrase XII / CA12 Protein (Catalog#50014-M08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse Carbonic Anhydrase XII / CA12 (rM Carbonic Anhydrase XII / CA12; Catalog#50014-M08H; NP_848483.2; Met 1-Ser 301).
Ig Type:	Rabbit IgG
Clone ID:	012
Specificity:	Mouse Carbonic Anhydrase XII / CA12
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA, IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:100-1:500
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse CA12.

Please Note: Optimal concentrations/dilutions should be determined by the end user.


Carbonic Anhydrase XII / Car12 Antibody, Rabbit MAb

Catalog Number: 50014-R012




Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Immunohistochemical staining of mouse CA12 in mouse kidney with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).



Immunohistochemical staining of mouse CA12 in mouse intestine with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).