

Cathepsin L/CTSL Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50015-T24

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse Cathepsin L/CTSL Protein (Catalog#50015-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse Cathepsin L/CTSL (rM Cathepsin L/CTSL; Catalog#50015-M08H; P06797; Met1-Asn334). Cathepsin L/CTSL specific IgG was purified by Mouse Cathepsin L/CTSL affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse Cathepsin L/CTSL
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA, IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse Cathepsin L/CTSL.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Cathepsin L/CTSL Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist

Catalog Number: 50015-T24



Immunohistochemical staining of mouse CTSL in mouse kidney with rabbit polyclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of mouse CTSL in mouse liver with rabbit polyclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.