



Catalog Number: 50020-R038

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse CD5L protein (Catalog#50020-M08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD5L / Aim (rM CD5L; Catalog#50020-M08H; Met 1-Val 352; NP_033820.2)
Ig Type:	Rabbit IgG
Clone ID:	38
Specificity:	Mouse CD5L / Aim
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse CD5L.

Please Note: Optimal concentrations/dilutions should be determined by the end user.