

CD86 Antibody (PE), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50068-R313-P

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CD86 Protein (Catalog#50068-M08H)
Reagents:	PE-conjugated Rabbit MAb antibody
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD86 (rM CD86; Catalog#50068-M08H; NP_062261.3; Met1-Glu245) and conjugated with PE under optimum conditions, the unreacted PE was removed.
Ig Type:	Rabbit IgG
Clone ID:	313
Specificity:	Mouse CD86
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>


CD86 Antibody (PE), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50068-R313-P

EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometric analysis of Mouse-CD86 expression on LPS-stimulated BABL/c splenocytes. Cells were stained with PE-conjugated anti-Mouse-CD86. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>