

CD132 / IL2RG Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50087-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse CD132 / IL2RG protein (Catalog#50087-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse CD132 / IL2RG (rM CD132 / IL2RG; Catalog#50087-M08H; NP_038591.1; Met 1-Ala 263). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse CD132 / IL2RG
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Mouse IL2RG.

Please Note: Optimal concentrations/dilutions should be determined by the end user.