

CD4 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50134-R766

EliteRmab® is a registered trademark of Sino Biological Inc.


GENERAL INFORMATION	
Immunogen:	Recombinant Mouse CD4 Protein (Catalog#50134-M08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD4 (rM CD4; Catalog#50134-M08H; NP_038516.1; Met1-Thr394).
Ig Type:	Rabbit IgG
Clone ID:	766
Specificity:	Mouse CD4
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:100-1:500

Please Note: Optimal concentrations/dilutions should be determined by the end user.




Catalog Number: 50134-R766


EliteRmab® is a registered trademark of Sino Biological Inc.



Immunohistochemical staining of mouse CD4 in mouse spleen with rabbit monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of mouse CD4 in mouse small intestine with rabbit monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of mouse CD4 in mouse colon with rabbit monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.