

ACVR2B / ACTRIIB Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50173-T16

GENERAL INFORMATION

Immunogen:	Recombinant Mouse ACVR2B protein (Catalog#50173-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse ACVR2B extracellular domain (rM ACVR2B; Catalog#50173-M08H; NP_031423.1; Met 1-Thr 134). ACVR2B specific IgG was purified by mouse ACVR2B affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse ACVR2B / ACTRIIB
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse ACVR2B.
-------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.