

Angiotensinogen/AGT Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50337-R001

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Mouse Angiotensinogen/AGT Protein (Catalog#50337-M08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse Angiotensinogen/AGT (rM Angiotensinogen/AGT; Catalog#50337-M08H; P11859; Met1-Val477).
Ig Type:	Rabbit IgG
Clone ID:	001
Specificity:	Mouse Angiotensinogen/AGT
Formulation:	0.2 µm filtered solution in PBS with 10% Trehalose, pH7.0
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB
----------------------	----

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:1000-1:5000
---------------------	-------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.


Angiotensinogen/AGT Antibody, Rabbit MAb

Catalog Number: 50337-R001



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Anti-AGT rabbit monoclonal antibody at
1:1000 dilution

Lane A: Mouse spleen tissue lysate

Lane B: Mouse liver tissue lysate

Lane C: Mouse kidney tissue lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000
dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 51 kDa