

ANGPTL4 / angiopoietin-like protein 4 Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50356-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse ANGPTL4 / angiopoietin-like protein 4 Protein (Catalog#50356-M07H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse ANGPTL4 / angiopoietin-like protein 4 (rM ANGPTL4 / angiopoietin-like protein 4; Catalog#50356-M07H; Q9Z1P8; Lys167-Ser410). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse ANGPTL4 / angiopoietin-like protein 4
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Mouse ANGPTL4 / angiopoietin-like protein 4.

Please Note: Optimal concentrations/dilutions should be determined by the end user.