

Factor IX/F9 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50362-T24

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse Factor IX/F9 Protein (Catalog#50362-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse Factor IX/F9 (rM Factor IX/F9; Catalog#50362-M08H; P16294; Met1-Thr471). Factor IX/F9 specific IgG was purified by Mouse Factor IX/F9 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse Factor IX/F9
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA, IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:200-1:1000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse Factor IX/F9.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Factor IX/F9 Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 50362-T24



Sino Biological
Biological Solution Specialist



Immunochemical staining of mouse F9 in mouse liver with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.



Immunochemical staining of mouse F9 in mouse kidney with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.