

CD8 alpha Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50389-T16

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CD8 alpha Protein (Catalog#50389-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse CD8 alpha (rM CD8 alpha; Catalog#50389-M08H; NP_001074579.1; Met1-Tyr196). CD8 alpha specific IgG was purified by Mouse CD8 alpha affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse CD8 alpha
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	BB154331, Ly-2, Ly-35, Ly-B, Lyt-2

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse CD8 alpha.
-------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.