

CD5 Antibody (PerCP), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50403-R249-C

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CD5 protein (Catalog#50403-M08H)
Reagents:	PerCP-conjugated Rabbit monoclonal antibody
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD5 (rM CD5; Catalog#50403-M08H; NP_031676.3; Met1-Pro371) and conjugated with PerCP under optimum conditions, the unreacted PerCP was removed.
Ig Type:	Rabbit IgG
Clone ID:	249
Specificity:	Mouse CD5
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Alternative Names:	CD5

APPLICATIONS

Applications:	FCM
----------------------	-----


RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.



Catalog Number: 50403-R249-C

EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometric analysis of mouse CD5 expression on splenocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.