

# Endoglin / CD105 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 50407-T24

| GENERAL INFORMATION       |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen:                | Recombinant Mouse Endoglin / CD105 Protein (Catalog#50407-M08H)                                                                                                                                                                                    |
| Preparation               | Produced in rabbits immunized with purified, recombinant Mouse Endoglin / CD105 (rM Endoglin / CD105; Catalog#50407-M08H; NP_031958.2; Met1-Gly580). Endoglin / CD105 specific IgG was purified by Mouse Endoglin / CD105 affinity chromatography. |
| Ig Type:                  | Rabbit IgG                                                                                                                                                                                                                                         |
| Specificity:              | Mouse Endoglin / CD105                                                                                                                                                                                                                             |
| Formulation:              | 0.2 µm filtered solution in PBS                                                                                                                                                                                                                    |
| Storage:                  | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles. |
| APPLICATIONS              |                                                                                                                                                                                                                                                    |
| Applications:             | ELISA, IHC-P                                                                                                                                                                                                                                       |
| RECOMMENDED CONCENTRATION |                                                                                                                                                                                                                                                    |
| IHC-P                     | IHC-P: 1:100-1:1000                                                                                                                                                                                                                                |
| ELISA                     | ELISA: 1:5000-1:10000<br>This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse Endoglin / CD105.                                                                                                     |

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**


# Endoglin / CD105 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 50407-T24

---



Immunohistochemical staining of mouse ENG in mouse liver with rabbit polyclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.