

CD55 / DAF Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50468-R084

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse CD55 / DAF protein (Catalog#50468-M08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD55 / DAF (rM CD55 / DAF; Catalog#50468-M08H; Q61475; Met1-Thr361).
Ig Type:	Rabbit IgG
Clone ID:	084
Specificity:	Mouse CD55 / DAF
	No cross-reactivity in ELISA with Human CD55
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA,ELISA(Cap)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse CD55 / DAF.
Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Mouse CD55/DAF in ELISA pair set (Catalog: # SEK50468). In a sandwich ELISA, it can be used as capture antibody.

Please Note: Optimal concentrations/dilutions should be determined by the end user.