

IFNAR-1 / CRF2-1 Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50469-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse IFNAR-1 protein (Catalog#50469-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse IFNAR-1 / CRF2-1 (rM IFNAR-1; Catalog#50469-M08H; NP_034638.2; Met 1-Thr 429). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse IFNAR-1 / CRF2-1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB, ELISA
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1:500-1:1000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse IFNAR-1.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

IFNAR-1 / CRF2-1 Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50469-RP01



Anti-IFNAR1 rabbit polyclonal antibody at 1:500 dilution
Lane A: K562 Membrane Lysate

Lysates/proteins at 20 µg per lane.
Secondary
Goat Anti-Rabbit IgG H&L (Dylight800) at 1/10000 dilution.

Developed using the Odyssey technique.
Performed under reducing conditions.

Predicted band size: 64 kDa
Observed band size: 55 kDa