

CD229 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50504-R016

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CD229 Protein (Catalog#50504-M08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD229 (rM CD229; Catalog#50504-M08H; NP_032560.2; Met1-Arg453).
Ig Type:	Rabbit IgG
Clone ID:	106
Specificity:	Mouse CD229
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA(Cap)
----------------------	------------

RECOMMENDED CONCENTRATION

ELISA(Cap)	ELISA(Cap): 1:250-1:2000 In a sandwich ELISA, the capture antibody (Catalog: # 50504-R016) can be used to detect Mouse CD229 when paired with the detection antibody (Catalog: # 50504-R001).
-------------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.