



Catalog Number: 50510-R004

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CD19 Protein (Catalog#50510-M08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD19 (rM CD19; Catalog#50510-M08H; NP_033974.2; Met1-Gly287).
Ig Type:	Rabbit IgG
Clone ID:	004
Specificity:	Human PDHA1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.


APPLICATIONS

Applications:	IHC-P
----------------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:100-1:500
--------------	--------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.



Immunohistochemical staining of Mouse CD19 in Mouse spleen with rabbit monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.