

# CD2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 50537-T16

## GENERAL INFORMATION

|                     |                                                                                                                                                                                                                             |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen:</b>   | Recombinant Mouse CD2 Protein (Catalog#50537-M08H)                                                                                                                                                                          |
| <b>Preparation</b>  | Produced in rabbits immunized with purified, recombinant Mouse CD2 (rM CD2; Catalog#50537-M08H; NP_038514.1; Met1-Ser203). CD2 specific IgG was purified by Mouse CD2 affinity chromatography.                              |
| <b>Ig Type:</b>     | Rabbit IgG                                                                                                                                                                                                                  |
| <b>Specificity:</b> | Mouse CD2                                                                                                                                                                                                                   |
| <b>Formulation:</b> | PBS, pH7.0 with 0.03% Proclin300                                                                                                                                                                                            |
| <b>Storage:</b>     | This antibody can be stored at 2℃-8℃ for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20℃ to -80℃. Avoid repeated freeze-thaw cycles. |

## APPLICATIONS

|                      |       |
|----------------------|-------|
| <b>Applications:</b> | ELISA |
|----------------------|-------|

## RECOMMENDED CONCENTRATION

|              |                                                                                                                                   |
|--------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>ELISA</b> | ELISA: 1:5000-1:10000<br>This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse CD2. |
|--------------|-----------------------------------------------------------------------------------------------------------------------------------|

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**