

Cathepsin E/CTSE Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50564-T24

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse Cathepsin E/CTSE Protein (Catalog#50564-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse Cathepsin E/CTSE (rM Cathepsin E/CTSE; Catalog#50564-M08H; NP_031825.2; Met1-Pro397). Cathepsin E/CTSE specific IgG was purified by Mouse Cathepsin E/CTSE affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse Cathepsin E/CTSE
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA,IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:200-1:1000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse Cathepsin E/CTSE.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Cathepsin E/CTSE Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist


Catalog Number: 50564-T24



Immunohistochemical staining of mouse CTSE in mouse liver with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of mouse CTSE in mouse stomach with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of mouse CTSE in mouse colon with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.