

Fetuin-B / FETUB Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50572-R110

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse Fetuin-B / FETUB protein (Catalog#50572-M08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse Fetuin-B / FETUB (rM Fetuin-B / FETUB; Catalog#50572-M08H; Q9QXC1-1; Met 1-Pro 388).
Ig Type:	Rabbit IgG
Clone ID:	110
Specificity:	Mouse Fetuin-B / FETUB
	No cross-reactivity in ELISA with Human FETUB
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Mouse FETUB.

Please Note: Optimal concentrations/dilutions should be determined by the end user.