

EphA2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50586-RP02

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse EphA2 protein (Catalog#50586-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse EphA2 (rM EphA2; Catalog#50586-M08H; NP_034269.2; Met 1-Asn 535). EphA2 specific IgG was purified by Mouse EphA2 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse EphA2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse EPHA2.

Please Note: Optimal concentrations/dilutions should be determined by the end user.