

E-Cadherin Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50671-T16

GENERAL INFORMATION

Immunogen:	Recombinant Mouse E-Cadherin Protein (Catalog#50671-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse E-Cadherin (rM E-Cadherin; Catalog#50671-M08H; P09803; Met1-Val709). E-Cadherin specific IgG was purified by Mouse E-Cadherin affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse E-Cadherin
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000
-------	-----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.