

CD59a Antibody (APC), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50724-R108-A

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CD59a Protein (Catalog#50724-M08H)
Reagents:	APC-conjugated Rabbit monoclonal antibody
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD59a (rM CD59a; Catalog#50724-M08H; NP_001104530.1; Met1-Lys95) and conjugated with APC under optimum conditions, the unreacted APC was removed.
Ig Type:	Rabbit IgG
Clone ID:	108
Specificity:	Mouse CD59a
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !
Alternative Names:	AA987121,Cd59,protectin,RP24-297H17.1

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>


CD59a Antibody (APC), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 50724-R108-A

EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometric analysis of Mouse CD59a expression on BABL/c bone marrow erythroid cells. Cells were stained with APC-conjugated anti-Mouse CD59a. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of erythroid cells.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>