

CD122/IL2RB Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50792-T16

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CD122/IL2RB Protein (Catalog#50792-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse CD122/IL2RB (rM CD122/IL2RB; Catalog#50792-M08H; P16297; Met1-Glu240). CD122/IL2RB specific IgG was purified by Mouse CD122/IL2RB affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse CD122/IL2RB
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000
-------	-----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.