

CNDP1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50959-T52

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CNDP1 Protein (Catalog#50959-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse CNDP1 (rM CNDP1; Catalog#50959-M08H; Q8BUG2; Met1-Tyr492). CNDP1 specific IgG was purified by Mouse CNDP1 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse CNDP1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, ELISA, IP
----------------------	---------------

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:2000
Immunoprecipitation	IP: 1-5 µL/mg of lysate
ELISA	ELISA: 1:5000-1:10000

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)


Website: <http://www.sinobiological.com>

CNDP1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50959-T52




Anti-CNDP1 rabbit polyclonal antibody at
1:500 dilution
Lane A: Mouse kidney tissue lysate

Lysates/proteins at 30 µg per lane.
Secondary
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000
dilution

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 57 kDa
Observed band size: 55 kDa



CNDP1 was immunoprecipitated using:
Lane A: 0.5 mg Mouse kidney Whole Cell
Lysate

4 µL anti-CNDP1 rabbit polyclonal antibody
and 60 µg of Immunomagnetic beads Protein
A/G.

Primary antibody:
Anti-CNDP1 rabbit polyclonal antibody, at
1:100 dilution

Secondary antibody:
Clean-Blot IP Detection Reagent (HRP) at
1:1000 dilution

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 57 kDa
Observed band size: 57 kDa

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>