

# IL12B Antibody, Mouse MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 51004-MM01

## GENERAL INFORMATION

|                     |                                                                                                                                                                                                                                                                                                                                       |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen:</b>   | Recombinant Mouse IL12B Protein (Catalog#51004-MCCH)                                                                                                                                                                                                                                                                                  |
| <b>Preparation</b>  | This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Mouse IL12B (rM IL12B; Catalog#51004-MCCH; P43432; Met1-Ser335). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography. |
| <b>Ig Type:</b>     | Mouse IgG1                                                                                                                                                                                                                                                                                                                            |
| <b>Clone ID:</b>    | 01                                                                                                                                                                                                                                                                                                                                    |
| <b>Specificity:</b> | Mouse IL12B                                                                                                                                                                                                                                                                                                                           |
| <b>Formulation:</b> | 0.2 µm filtered solution in PBS                                                                                                                                                                                                                                                                                                       |
| <b>Storage:</b>     | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.                                                                                    |


## APPLICATIONS

|                      |       |
|----------------------|-------|
| <b>Applications:</b> | IHC-P |
|----------------------|-------|

## RECOMMENDED CONCENTRATION

|              |                    |
|--------------|--------------------|
| <b>IHC-P</b> | IHC-P: 1:100-1:500 |
|--------------|--------------------|

**Please Note:** Optimal concentrations/dilutions should be determined by the end user.



Immunohistochemical staining of Mouse IL12B in Mouse spleen with rabbit polyclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.

Manufactured By Sino Biological Inc., FOR RESEARCH USE ONLY. NOT FOR USE IN HUMANS.

Fax :+86-10-58628220 ● Tel:+86-400-890-9989 ● <http://www.sinobiological.com>