

IL-12 p40/IL12B Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 51004-R240

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse IL-12 p40/IL12B Protein (Catalog#51004-MCCH)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse IL-12 p40/IL12B (rM IL-12 p40/IL12B; Catalog#51004-MCCH; P43432; Met1-Ser335).
Ig Type:	Rabbit IgG
Clone ID:	240
Specificity:	Mouse IL-12 p40/IL12B
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000

Please Note: Optimal concentrations/dilutions should be determined by the end user.