



Catalog Number: 51048-R017

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse GM-CSF/CSF2 Protein (Catalog#51048-M07H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse GM-CSF/CSF2 (rM GM-CSF/CSF2; Catalog#51048-M07H; P01587; Ala18-Lys141).
Ig Type:	Rabbit IgG
Clone ID:	017
Specificity:	Mouse GM-CSF/CSF2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA, ELISA(Cap)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse GM-CSF/CSF2.
Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Mouse GM-CSF/CSF2 in ELISA pair set (Catalog: # SEK51048). In a sandwich ELISA, it can be used as capture antibody.

Please Note: Optimal concentrations/dilutions should be determined by the end user.