

HER4 / ErbB4 Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 51064-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse ErbB4 / HER4 protein (Catalog#51064-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse ErbB4 / HER4 (rM ErbB4 / HER4; Catalog#51064-M08H; Q99P91; Met1-Leu652). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse ErbB4 / HER4
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Mouse ErbB4 / HER4.

Please Note: Optimal concentrations/dilutions should be determined by the end user.