

IL-17A Antibody (FITC), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 51065-R208-F

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Mouse IL-17A Protein (Catalog#51065-MNAE)
Reagents:	FITC-conjugated Rabbit MAb antibody
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse IL-17A (rM IL-17A; Catalog#51065-MNAE; Q62386; Thr22-Ala158) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.
Ig Type:	Rabbit IgG
Clone ID:	208
Specificity:	Mouse IL-17A
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>


IL-17A Antibody (FITC), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 51065-R208-F

EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometric analysis of Mouse IL17 expression on EL4 cells. The cells were stimulated for 4-6 hours with PMA and Ionomycin in the presence of GolgiPlug. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), stained with FITC-conjugated anti-Mouse IL17. The dot plots were derived from gated of intact cells.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>