

IL36B/IL1F8 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 51071-T24

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse IL36B/IL1F8 Protein (Catalog#51071-MNAE)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse IL36B/IL1F8 (rM IL36B/IL1F8; Catalog#51071-MNAE; Q9D6Z6; Met1-Lys183). IL36B/IL1F8 specific IgG was purified by Mouse IL36B/IL1F8 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse IL36B/IL1F8
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA,IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse IL36B/IL1F8.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

IL36B/IL1F8 Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 51071-T24



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of mouse IL1F8 in mouse small intestine with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).



Immunohistochemical staining of mouse IL1F8 in mouse stomach with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).