

CD8 / CD8 alpha / Leu-2 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 60001-R070

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Ferret CD8 / CD8 alpha / Leu-2 Protein (Catalog#60001-F08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Ferret CD8 / CD8 alpha / Leu-2 (Catalog#60001-F08H; ABS50091.1; Met1-Glu186).
Ig Type:	Rabbit IgG
Clone ID:	070
Specificity:	Ferret CD8 / CD8 alpha / Leu-2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA(Cap)
RECOMMENDED CONCENTRATION	
Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Ferret CD8 / CD8 alpha / Leu-2 in ELISA pair set (Catalog: # SEK60001). In a sandwich ELISA, it can be used as capture antibody when paired with (Catalog: # 60001-MM02).

Please Note: Optimal concentrations/dilutions should be determined by the end user.