

CD45RA Antibody (PerCP), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 68064-MM27-C

GENERAL INFORMATION

Reagents:	PerCP-conjugated Mouse monoclonal antibody
Ig Type:	Mouse IgG1
Clone ID:	27
Specificity:	Human CD45RA
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2℃-8℃ for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !
Alternative Names:	CD45RA

APPLICATIONS

Applications:	FCM
---------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD45RA Antibody (PerCP), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 68064-MM27-C



Flow cytometric analysis of Human CD45RA expression on human whole blood lymphocytes. Cells were stained with PerCP-conjugated anti-Human CD45RA. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.