

Anti-Integrin $\alpha 4\beta 7$ (Research Grade Vedolizumab Biosimilar)



Sino Biological
Biological Solution Specialist

Catalog Number: 68121-H001


GENERAL INFORMATION

Preparation	This product is a recombinant monoclonal antibody expressed from HEK293 cells.
Ig Type:	IgG1, Kappa
Formulation:	Histidine and Arginine buffer containing 120mM NaCl, 0.02% Tween 80, pH6.0
Purity:	≥ 95 % as determined by SDS-PAGE. ≥ 95 % as determined by SEC-HPLC.
Endotoxin:	< 1 EU/mg
Storage:	Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C . Preservative-Free. Avoid repeated freeze-thaw cycles.


RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.

SDS-PAGE



SEC-HPLC



For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>