

Anti-FZD1/2/5/7/8 (Research Grade Vantictumab Biosimilar)



Sino Biological
Biological Solution Specialist

Catalog Number: 68204-H001


GENERAL INFORMATION

Preparation	This product is a recombinant monoclonal antibody expressed from HEK293 cells.
Ig Type:	IgG2, Lambda
Formulation:	Histidine and Arginine buffer containing 120mM NaCl, 0.02% Tween 80, pH6.0
Purity:	≥95 % as determined by SDS-PAGE. ≥95 % as determined by SEC-HPLC.
Endotoxin:	< 1 EU/mg
Storage:	Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.


RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.

SDS-PAGE



SEC-HPLC



For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>