

# Anti-HMPV-F (Research Grade Clone MPV458)



Sino Biological  
Biological Solution Specialist

Catalog Number: 68251-H001A


## GENERAL INFORMATION

|                     |                                                                                                                                                          |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Preparation</b>  | This product is a recombinant monoclonal antibody expressed from CHO cells.                                                                              |
| <b>Ig Type:</b>     | IgG1, Kappa                                                                                                                                              |
| <b>Clone ID:</b>    | MPV458                                                                                                                                                   |
| <b>Formulation:</b> | Histidine and Arginine buffer containing 120mM NaCl, 0.02% Tween 80, pH6.0                                                                               |
| <b>Purity:</b>      | ≥95 % as determined by SDS-PAGE.<br>≥95 % as determined by SEC-HPLC.                                                                                     |
| <b>Endotoxin:</b>   | < 1 EU/mg                                                                                                                                                |
| <b>Storage:</b>     | Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles. |


## RECOMMENDED CONCENTRATION

*Please Note: Optimal concentrations/dilutions should be determined by the end user.*

### SDS-PAGE



### SEC-HPLC



**For Research Use Only. Not for use in diagnostic or therapeutic procedures.**

**Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)**

**Website: <http://www.sinobiological.com>**