

Anti-DM1 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 68254-R0002

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	DM1 conjugated to BSA
Preparation:	This product is an antibody recombinantly expressed in HEK293 cells.
Ig Type:	Rabbit IgG
Clone ID:	0002
Specificity:	Specifically recognize DM1.
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 0.01-0.1 µg/mL
--------------	-----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.


Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>




Catalog Number: 68254-R0002

EliteRmab® is a registered trademark of Sino Biological Inc.



Immobilized DM1-OVA at 1 µg/mL (100 µL/well), free DM1 or DM4 and Anti-DM1 Antibody add to the plate (total volume 50 µL/well, Anti-DM1 Antibody working concentration is 0.005 µg/mL), the binding of DM1-OVA and Anti-DM1 Antibody were blocked by the gradient-diluted free DM1 or DM4 from 0.07 ng/mL.



Immobilized DM1-OVA at 0.1 µg/mL (100 µL/well) can bind the Anti-DM1 Antibody, the EC₅₀ is 8.53 ng/mL.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>