

Anti-PD1 (Research Grade Nofazinlimab Biosimilar)



Sino Biological
Biological Solution Specialist

Catalog Number: 68279-H001A

GENERAL INFORMATION

Preparation	This product is a recombinant monoclonal antibody expressed from HEK293 cells.
Ig Type:	IgG1, Kappa
Formulation:	Histidine and Arginine buffer containing 120mM NaCl, 0.02% Tween 80, pH6.0
Purity:	≥95 % as determined by SDS-PAGE. ≥95 % as determined by SEC-HPLC.
Endotoxin:	< 1 EU/mg
Storage:	Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.